

Latest BS3+

Wi-Fi Digital Microscope Dermatoscope



www.portable-ultrasound-scanner


Quality, Professional, Innovational

Product Specifications

Pixels	2,000,000 pixels
Enlarge multiple	50~1000X
Video pixels	1920*1080P/1280*720P/640*480P
Imaging distance	Manual adjustment (3~60mm)
Wi-Fi connection	10m (open environment)
Image format	JPG
USB interface type	USB2.0
Power supply	USB or 5V DC
Video format	MP4/AVI
Dynamic frames	15f/s under 600LUX Brightness
Operation system	Android 4.3 , iOS 8.0
Size	142mm(L)*37mm(R)
Net weight	80g
Standard parts list	Wi-Fi microscope User manual 1m data cable Plastic base Metal tripod Packing box



Images



Actual product pictures



Wuxi Biomedical Technology Co., Ltd. www.portable-ultrasound-scanner.com

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbiomedical.com

Mobile phone: 86-13382885910

BIO is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.